

## Standard Operating Procedure:

# Gold Etch

### Table of Contents:

<b>Table of Contents:</b> .....	<b>1</b>
<b>Principle of Operation</b> .....	<b>2</b>
<b>Material Requirements</b> .....	<b>2</b>
Equipment:.....	2
Chemicals:.....	2
Personal Protective Equipment:.....	2
<b>Procedure</b> .....	<b>2</b>
Gold Etch.....	2
DI Water Rinse.....	3
Sample Dry.....	3
Cleanup.....	3
<b>Accident Procedure</b> .....	<b>3</b>
Contact.....	3
Spills.....	4
<b>Revision History:</b> .....	<b>4</b>

## Principle of Operation

To remove gold from a substrate – patterned or unpatterned.

## Material Requirements

### Equipment:

Substrate, glass Petri dish or beaker (for etchant), glass beaker (for rinse) and stainless steel tweezers

### Chemicals:

Gold Etchant (contains Iodine, Potassium Iodide and water)

- Gold Etchant Hazards
  - o Vapors severely irritate the eyes, mucous membranes and respiratory tract.
  - o Do not concentrate the etchant by evaporating off the water. Iodine and Potassium Iodide in concentrated forms can cause severe skin irritation in concentrated solutions or crystalline form.
  - o Do not mix with strong reducing agents, ammonia, ammonium hydroxide, powdered metals or alkali metals. Toxic by products include oxides of iodine and iodine fumes.

### Personal Protective Equipment:

Trionic gloves on top of nitrile gloves, apron, safety glasses and face-shield

## Procedure

### Gold Etch

1. Rinse both beakers with DI water prior to beginning the process.
2. Stand the beaker to be used for rinsing on a few fab wipes in the hood and fill it with DI water so that the water level will cover the entire substrate.
3. Get a glass beaker that will fit your samples for processing – you should find one labeled “Gold Etchant” on the shelves. Place it on fab wipes in the hood.
4. Carefully pour some of the Gold Etchant into the glass beaker so that there is enough liquid to cover your submerged sample.
5. Calculate the etch time for your sample. You will need to know the thickness of your gold layer. At 25°C, the gold will etch at a rate of approximately 28 Å/sec.
6. Put your sample into the etchant and soak for the appropriate amount of time calculated in the previous step. Agitation of the Petri dish will accelerate the etch.

## DI Water Rinse

1. When the etch is complete, transfer the sample carefully to the DI water rinse beaker.
2. If you used tweezers to move the sample, make sure you leave them in the rinse beaker as well.
3. Let the sample and tools soak in DI water for 5 minutes.
4. Remove the sample from the rinse container.
5. Rinse the sample with DI water in the hood sink.

## Sample Dry

1. After the water rinse is finished, blow the sample dry with the N<sub>2</sub> gun.
2. After getting most of the water off, you can dry the samples more in an oven or on a hotplate if allowable for your sample.
3. Inspect sample for traces of un-etched gold. If features are small, use an optical microscope. If more etch time is required, place wafer back into Petri dish with the etchant for another 30 seconds while agitating. Repeat rinse and drying procedure.

## Cleanup

1. The etchant may be used for multiple etches. For temporary storage (<1 day), place the top of the Petri dish over the etchant and store on fab wipes in the back of the hood. Make sure the dish is clearly labeled "Gold Etchant".
2. Dump the etchant waste into the carboy designated for acids.
3. Rinse the Petri dish once with DI water, and dump it into the same carboy.
4. Dump the DI rinse beakers into the acid carboy.
5. Rinse all the containers again with DI water in the hood sink. This time, dump them into the sink within the process hood.
6. Return all labware to its proper location. The Petri dish and the beaker can drip dry on fab wipes in the hood; however, remember to move them back to their storage location once dry.
7. Wipe up any drips in the area with chemical wipes and dispose in the trash.
8. Store the Gold Etchant in the oxidizers cabinet.
9. Inspect all of the PPE to ensure it did not come in contact with the etchant before returning it to its storage location.

## Accident Procedure

### Contact

- Skin: Remove contaminated clothing, rinse affected area with water for at least 15 minutes. Wash skin with soap and water. **If there is any irritation, get immediate medical attention.**
- Eye: Immediately flush with water for at least 15 minutes while lifting upper and lower eyelids occasionally. **Get immediate medical attention.**

- Ingestion: Drink water or milk unless victim is unconscious. Induce vomiting if help is not immediately available. **Get immediate medical attention.**
- Inhalation: Remove to fresh air. Resuscitate if necessary. Take care not to inhale any fumes released from the victim's lungs. **Get immediate medical attention.**

## Spills

If a small, contained spill occurs, such as inside the hood, wipe it up with chemical wipes and dispose of them in the proper trash container. If a large spill occurs, evacuate the area and notify the cleanroom staff.

## Revision History:

- Version 0.0 – Created September 8<sup>th</sup>, 2017
- Version 1.0 – Revised May 22, 2025